

Tasso+ Collection Kits



Click. Collect. Done!

Blood collection has never been this easy. Tasso+ Kits contain everything a patient needs to collect their own blood sample.




Create a Better Patient Experience

Collecting samples with Tasso+ is convenient and virtually painless, and eliminates common barriers to completion.



Ship Samples with Confidence


Tasso+ Kits are designed to ship samples from the patient to the lab, enabling remote and hybrid workflows with ease.




Accurately Measure Blood Analytes

Tasso+ Kits enable the collection of high volume, clinical grade blood samples with multiple standard preparation options.


Collection Process



1 Receive kit containing all supplies for sample collection.



2 Collect patient sample.



3 Place sample in provided packaging.



4 Return sample for analysis.

Technical Specifications

Indication for Use: Allows self-collection of a capillary blood sample in conjunction with a compatible sample collection reservoir that is shipped to an accredited laboratory for analysis.

Kit Dimensions: 180x157x38mm (7.1x6.2x1.5in)

Kit Weight: 178g (6.25oz)

Shipping: Sample can be shipped using packaging materials provided with the kit.

Kit Storage: 18-25° C (64.4 - 77° F)

Shelf Life: Up to 18 months, dependent on kit components. Tasso+ device is 24 months.

Availability: Kits are procedure packs per Article 22 of the MDR (EU) 2017/745. The Tasso+ device is CE Marked and MHRA registered.



Safety Specifications

Biocompatibility: ANSI AAMI ISO 10993

Sterilization: ANSI AAMI ISO 11137
ANSI AAMI ISO 11607

Kit Configurations

Product	Base Part Number	Tube Type	Blood Collection Vol.	Blood Collection Additive
Tasso+ EDTA Kit	K1025-01-EDT-G02-APP	K ₂ EDTA (lavender)	250-500µL	K ₂ EDTA
Tasso+ Serum Gel Kit	K1025-01-SGT-G02-APP	Serum Separator (gold)	400-600µL	Clot Activator with Gel
Tasso+ LH Kit	K1025-01-LHT-G02-APP	Lithium Heparin (green)	200-400µL	Lithium Heparin
Tasso+ Plasma Gel Kit	K1025-01-PGT-G02-APP	Plasma Separator (light green)	400-600µL	Lithium Heparin with Gel

Please visit www.tassoinc.com/eu/tasso-plus-kit for more information.