


# Tasso+ Collection Kits



Research Use Only

## Click. Collect. Done!

Blood collection has never been this easy. Tasso+ Kits contain everything a patient needs to collect their own blood sample.




### Create a Better Patient Experience

Collecting samples with Tasso+ is convenient and virtually painless, and eliminates common barriers to completion.



### Ship Samples with Confidence


Tasso+ Kits are designed to ship samples from the patient to the lab, enabling remote and hybrid workflows with ease.




### Accurately Measure Blood Analytes

Tasso+ Kits enable the collection of high volume, clinical grade blood samples with multiple standard preparation options.


## Collection Process



1 Receive kit containing all supplies for sample collection.



2 Collect patient sample.



3 Place sample in provided packaging.



4 Return sample for analysis.

## Technical Specifications

**Indication for Use:** Allows self-collection of microliter capillary whole blood samples in a compatible sample collection reservoir that is shipped to an accredited laboratory for analysis. For Research Use Only. Not for use in diagnostic procedures.

**Kit Dimensions:** 7.1x6.2x1.5in

**Kit Weight:** 6.25oz/178g

**Shipping:** Sample can be shipped using packaging materials provided with the kit.

**Kit Storage:** 18-25° C (64.4 - 77° F)

**Shelf Life:** Up to 18 months, dependent on kit components. Tasso+ device is 24 months.

**Availability:** Kits are Research Use Only in the United States. Tasso+ device is Class II Lancet 510(k)-Cleared.



## Safety Specifications

**Biocompatibility:** ANSI AAMI ISO 10993

**Sterilization:** ANSI AAMI ISO 11137  
ANSI AAMI ISO 11607

## Kit Configurations

| Product               | Base Part Number     | Tube Type                      | Blood Collection Vol. | Blood Collection Additive |
|-----------------------|----------------------|--------------------------------|-----------------------|---------------------------|
| Tasso+ EDTA Kit       | K1025-01-EDT-G01-RUO | K <sub>2</sub> EDTA (lavender) | 250-500µL             | K <sub>2</sub> EDTA       |
| Tasso+ Serum Gel Kit  | K1025-01-SGT-G01-RUO | Serum Separator (gold)         | 400-600µL             | Clot Activator with Gel   |
| Tasso+ LH Kit         | K1025-01-LHT-G01-RUO | Lithium Heparin (green)        | 200-400µL             | Lithium Heparin           |
| Tasso+ Plasma Gel Kit | K1025-01-PGT-G01-RUO | Plasma Separator (light green) | 400-600µL             | Lithium Heparin with Gel  |

Please visit [www.tassoinc.com/tasso-plus-kit](http://www.tassoinc.com/tasso-plus-kit) for more information.